

## Product Datasheet

# Methyl-HP1 gamma/CBX3-H175 antibody (orb1149895)

<b>Description</b>	Rabbit polyclonal antibody to Methyl-HP1 gamma/CBX3-H175
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC
<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human HP1 gamma/CBX3(NP_009207.2).
<b>Preservatives</b>	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	20kDa
<b>Uniprot ID</b>	<a href="#">Q13185</a>
<b>Dilution Range</b>	WB 1:500 - 1:2000
<b>Expiration Date</b>	12 months from date of receipt.