

## Product Datasheet

# NMDAR1 (Phospho-S896/S897) antibody (orb1149957)

<b>Description</b>	Rabbit polyclonal antibody to NMDAR1 (Phospho-S896/S897)
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A synthetic phosphorylated peptide around S896 & S897 of human NMDAR1(NP_015566.1).
<b>Preservatives</b>	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	99kDa/101kDa/103kDa/105kDa/106kDa/107kDa
<b>Uniprot ID</b>	<a href="#">Q05586</a>
<b>Dilution Range</b>	WB 1:500 - 1:2000
<b>Expiration Date</b>	12 months from date of receipt.