

Product Datasheet

FGFR1 (Phospho-Y653/654) antibody (orb1149995)

Description	Rabbit polyclonal antibody to FGFR1 (Phospho-Y653/654)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y653 & Y654 of human FGFR1 (NP_075598.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	150kDa
Uniprot ID	P11362
Entrez	2260
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.