

Product Datasheet

SARS-CoV-2 NSP8 antibody (orb1150033)

Description	Rabbit polyclonal antibody to SARS-CoV-2 NSP8
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of coronavirus NSP8(YP_009742615.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25KDa
Uniprot ID	P0DTC2
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.