

Product Datasheet

HCoV-HKU1 Spike S1 antibody (orb1150057)

Description	Rabbit polyclonal antibody to HCoV-HKU1 Spike S1
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of coronavirus Spike S1(YP_173238.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75KDa/110KDa/145KDa
Uniprot ID	Q5MQD0
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.