

Product Datasheet

HCoV-OC43 Spike S2 antibody (orb1150220)

Description	Rabbit polyclonal antibody to HCoV-OC43 Spike S2
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 782-881 of coronavirus Spike S2 (YP_009555241.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75kDa/145kDa
Uniprot ID	P36334
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.