

Product Datasheet

MAPK14 (Phospho-T180/Y182) +MAPK11 (Phospho-T180/Y182) +MAPK13 (Phospho-T180/Y182) +MAPK12 (Phospho-T183/Y1845) antibody (orb1150252)

Description	Rabbit polyclonal antibody to MAPK14 (Phospho-T180/Y182) +MAPK11 (Phospho-T180/Y182) +MAPK13 (Phospho-T180/Y182) +MAPK12 (Phospho-T183/Y1845)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T180 & Y182 of human MAPK14MAPK11(NP_002742.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43KDa
Uniprot ID	Q15759 , Q16539 , O15264 , P53778
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.