

Product Datasheet

PDGFRB (Phospho-Y1009) antibody (orb1150259)

Description	Rabbit polyclonal antibody to PDGFRB (Phospho-Y1009)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1009 of human PDGFRB(NP_002600.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	190kDa
Uniprot ID	P09619
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.