

Product Datasheet

MEK1 (Phospho-T292) antibody (orb1150302)

Description	Rabbit polyclonal antibody to MEK1 (Phospho-T292)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T292 of human MEK1(NP_002746.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45KDa
Uniprot ID	Q02750
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.