

## Product Datasheet

### HCoV-HKU1 Spike S2 antibody (orb1150380)

<b>Description</b>	Rabbit polyclonal antibody to HCoV-HKU1 Spike S2
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	A synthetic peptide corresponding to a sequence within amino acids 800-900 of coronavirus Spike S2 (YP_173238.1).
<b>Preservatives</b>	PBS with 0.05% proclin300,50% glycerol,pH7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	75kDa/145kDa
<b>Uniprot ID</b>	<b>Q5MQD0</b>
<b>Dilution Range</b>	WB 1:500 - 1:2000
<b>Expiration Date</b>	12 months from date of receipt.