

Product Datasheet

SARS-CoV-2 NSP10 antibody (orb1150397)

Description	Rabbit polyclonal antibody to SARS-CoV-2 NSP10
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of coronavirus NSP10(YP_009725306.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20-25KDa
Uniprot ID	P0DTD1
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.