

Product Datasheet

IRE1 (Phospho-Ser724/726) antibody (orb1150426)

Description	Rabbit polyclonal antibody to IRE1 (Phospho-Ser724/726)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around of human Phospho-IRE1-Ser724/726.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O75460
Entrez	2081
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.