

Product Datasheet

Methyl-hnRNP A1-H173 antibody (orb1150477)

Description	Rabbit polyclonal antibody to Methyl-hnRNP A1-H173
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human Methyl-hnRNP A1-H173, 3-methyl (tau-methyl) .
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P09651
Entrez	3178
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.