

Product Datasheet

CDAN1 antibody [Out of stock] (orb1150601)

Description	Rabbit polyclonal antibody to CDAN1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human CDAN1(NP_612486.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IWY9
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.