

Product Datasheet

SARS-CoV-2 Spike S1 antibody (orb1150669)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Spike S1
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	DOT, ICC, IF, IP, WB
Immunogen	Recombinant fusion protein coronavirus Spike S1(YP_009724390.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110KDa
Uniprot ID	P0DTC2
Dilution Range	WB 1:500 - 1:2000 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.