

Product Datasheet

p53 (Phospho-T81) antibody (orb1150676)

Description	Rabbit polyclonal antibody to p53 (Phospho-T81)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T81 of human TP53
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	53KDa
Uniprot ID	P04637
Entrez	7157
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.