

## Product Datasheet

### MLKL (Phospho-S345) antibody (orb1150706)

<b>Description</b>	Rabbit polyclonal antibody to MLKL (Phospho-S345)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A phospho synthetic peptide of Mouse MLKL.
<b>Preservatives</b>	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>MW</b>	54KDa
<b>Uniprot ID</b>	<a href="#">Q9D2Y4</a>
<b>Entrez</b>	<a href="#">74568</a>
<b>Dilution Range</b>	WB 1:500 - 1:2000
<b>Expiration Date</b>	12 months from date of receipt.