

Product Datasheet

SARS-CoV-2 NSP13 antibody (orb1150786)

Description	Rabbit polyclonal antibody to SARS-CoV-2 NSP13
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of SARS-CoV-2 NSP13.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	75KDa
Uniprot ID	P0DTD1
Entrez	43740578
Dilution Range	WB 1:500 - 1:2000 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.