

Product Datasheet

Histone H3 (Phospho-T45) antibody (orb1150792)

Description	Rabbit polyclonal antibody to Histone H3 (Phospho-T45)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T45 of human Histone H3.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	17kDa
Uniprot ID	P68431
Entrez	8350
Dilution Range	WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.