

## Product Datasheet

### TrkB (Phospho-Y705) antibody (orb1150807)

<b>Description</b>	Rabbit polyclonal antibody to TrkB (Phospho-Y705)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A phospho specific peptide corresponding to residues surrounding Y705 of human Phospho-TrkB-Y705
<b>Preservatives</b>	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>MW</b>	140kDa
<b>Uniprot ID</b>	<a href="#">Q16620</a>
<b>Entrez</b>	<a href="#">4915</a>
<b>Dilution Range</b>	WB 1:500 - 1:2000
<b>Expiration Date</b>	12 months from date of receipt.