

Product Datasheet

ULK1 (Phospho-S467) antibody (orb1150834)

Description	Rabbit polyclonal antibody to ULK1 (Phospho-S467)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S467 of human ULK1
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	140KDa
Uniprot ID	O75385
Entrez	8408
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.