

Product Datasheet

5-HT-TriMethyl-Histone H3-K4 antibody (orb1150900)

Description	Rabbit polyclonal antibody to 5-HT-TriMethyl-Histone H3-K4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human 5-HT-TriMethyl-Histone H3-K4.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P68431
Entrez	8350
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.