

Product Datasheet

IKK β (Phospho-Y199) antibody (orb1151005)

Description	Rabbit polyclonal antibody to IKK β (Phospho-Y199)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y199 of human IKK β (NP_001547.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	87KDa
Uniprot ID	O14920
Entrez	3551
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.