

Product Datasheet

p70 S6 Kinase 1 (Phospho-T421) antibody (orb1151021)

Description	Rabbit polyclonal antibody to p70 S6 Kinase 1 (Phospho-T421)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T421 of human P70S6K1
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P23443
Entrez	6198
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.