



Product Datasheet

GM130 antibody (orb1151045)


Description	Rabbit polyclonal antibody to GM130
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human GM130(NP_004477.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130KDa
Uniprot ID	Q08379
Dilution Range	IHC-P: 1:50 - 1:200, IF/ICC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of PC-12 cel...