

Product Datasheet

SARS-CoV-2 Nucleoprotein antibody (orb1151059)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Nucleoprotein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	DOT, ELISA, IF, IP, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of coronavirus Nucleoprotein (P0DTC9).
Preservatives	PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Dilution Range	WB 1:500 - 1:2000 IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.