

Product Datasheet

AMPK β 1 (Phospho-S182) antibody (orb1151157)

Description	Rabbit polyclonal antibody to AMPK β 1 (Phospho-S182)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S182 of human AMPK β 1(NP_006244.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38KDa
Uniprot ID	Q9Y478
Dilution Range	WB 1:500 - 1:2000 IHC 1:50 - 1:200
Expiration Date	12 months from date of receipt.