

Product Datasheet

DiMethyl-Histone H3-K79 antibody (orb1151166)

Description	Rabbit polyclonal antibody to DiMethyl-Histone H3-K79
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	A synthetic dimethylated peptide around K79 of human histone H3(NP_003520.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17KDa
Uniprot ID	Q16695
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.