

Product Datasheet

IKK α (Phospho-T23) antibody (orb1151171)

Description	Rabbit polyclonal antibody to IKK α (Phospho-T23)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T23 of human IKK α .
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O15111
Entrez	1147
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.