

Product Datasheet

SARS-CoV-2 Spike S2 antibody (orb1151188)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Spike S2
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of SARS-CoV-2 Spike S2(YP_009724390.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36KDa
Uniprot ID	P0DTC2
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.