

Product Datasheet

SARS-CoV-2 M₁ protein antibody (orb1151198)

Description	Rabbit polyclonal antibody to SARS-CoV-2 M ₁ protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of SARS-CoV-2 M ₁ protein.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P0DTC5
Entrez	43740571
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.