

Product Datasheet

HCoV-NL63 Spike S2 antibody (orb1151213)

Description	Rabbit polyclonal antibody to HCoV-NL63 Spike S2
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1200-1300 of coronavirus Spike S2 (YP_003767.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160kDa
Uniprot ID	Q6Q1S2
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.