

Product Datasheet

KMT2D antibody (orb1151250)

Description	Rabbit polyclonal antibody to KMT2D
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human KMT2D(NP_003473.3).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14686
Dilution Range	IHC 1:50 - 1:200
Expiration Date	12 months from date of receipt.