

Product Datasheet

MERS-CoV Spike RBD antibody (orb1151282)

Description	Rabbit polyclonal antibody to MERS-CoV Spike RBD
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of coronavirus Spike RBD(YP_009047204.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	200KDa
Uniprot ID	K9N5Q8
Dilution Range	WB 1:500 - 1:2000IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.