

Product Datasheet

ErbB3/HER3 (Phospho-Y1289) antibody (orb1151328)

Description	Rabbit polyclonal antibody to ErbB3/HER3 (Phospho-Y1289)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1289 of human ERBB3 (NP_001973.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P21860
Entrez	2065
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.