

Product Datasheet

SARS-CoV-2 Spike antibody (orb1151340)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to SARS-CoV-2 Spike |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | DOT, ICC, IF, IP, WB |
| Immunogen | A synthetic peptide of coronavirus Spike(YP_009724390.1). |
| Preservatives | PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 110KDa/180KDa |
| Uniprot ID | PODTC2 |
| Dilution Range | WB 1:500 - 1:2000IF 1:50 - 1:200 |
| Expiration Date | 12 months from date of receipt. |