


Product Datasheet

GATA3 antibody (orb1151475)


Description	Rabbit monoclonal antibody to GATA3
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A conserved peptide between the GATA trans-activation and DNA-binding domain was used as the immunogen for the recombinant GATA3 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the recombinant GATA3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant GATA3 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	P23771
Hazard Information	This recombinant GATA3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml, Western blot: 1-2ug/ml
Expiration Date	12 months from date of receipt.



Western blot testing of human MCF7 cell ...



IHC staining of FFPE human bladder carci...



Negative control: IHC staining of FFPE h...