



Product Datasheet

NeuN antibody (orb1151476)

Description	Rabbit monoclonal antibody to NeuN
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human RNA binding protein fox-1 homolog 3 protein was used as the immunogen for the recombinant NeuN antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the recombinant NeuN antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant NeuN antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	A6NFN3
Hazard Information	This recombinant NeuN antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml, Western blot: 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human cerebellum ti...



IHC staining of FFPE human cerebrum tiss...



Western blot testing of human brain tiss...