



Product Datasheet

p40 antibody (orb1151485)


Description	Rabbit monoclonal antibody to p40
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IHC-P
Immunogen	A portion of amino acids 1-100 was used as the immunogen for the recombinant deltaNp63 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the recombinant deltaNp63 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant deltaNp63 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Purity	Protein A affinity
Uniprot ID	Q9H3D4
Hazard Information	This recombinant deltaNp63 antibody is available for research use only.
Dilution Range	ELISA (order BSA-free format for coating), Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate tiss...



IHC staining of FFPE human prostate tiss...



FACS staining of PFA-fixed human HeLa ce...