

Product Datasheet

CD48 (Pan Leukocyte Marker) antibody (orb1151520)

Description

Mouse monoclonal antibody to CD48 (Pan Leukocyte



Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

A portion of amino acids 59-178 was used as the immunogen for the CD48 antibody.

IHC staining
of FFPE
human
tonsil
tissue...

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Aliquot the CD48 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note


For research use only

Application notes

Optimal dilution of the CD48 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide



IHC staining
of FFPE
human
lymph node
ti...

Isotype

Mouse IgG2c, kappa

Clonality

Monoclonal

Purity

Protein A/G affinity

Uniprot ID

P09326


Hazard

Information

This CD48 antibody is available for research use only.

Dilution Range

Western blot: 2-4ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml



Western blot
testing of
human Raji
cell ...

Expiration Date

12 months from date of receipt.