



Product Datasheet

CD5 antibody (orb1151521)


Description	Mouse monoclonal antibody to CD5
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length protein corresponding to human CD5 was used as the immunogen for the recombinant CD5 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the recombinant CD5 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant CD5 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	P06127
Hazard Information	This recombinant CD5 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml, Western blot: 1-2ug/ml
Expiration Date	12 months from date of receipt.



Western blot testing of human Jurkat cel...



IHC staining of FFPE human tonsil tissue...



Negative control: IHC staining of FFPE h...