



## Product Datasheet

### Prolactin antibody (orb1151538)


|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Prolactin                                                                       |
| <b>Species/Host</b>        | Mouse                                                                                                        |
| <b>Reactivity</b>          | Human                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                 |
| <b>Tested Applications</b> | IHC-P                                                                                                        |
| <b>Immunogen</b>           | A portion of amino acids 63-201 was used as the immunogen for the recombinant PRL antibody.                  |
| <b>Preservatives</b>       | 0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide                                 |
| <b>Storage</b>             | Aliquot the recombinant PRL antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                        |
| <b>Application notes</b>   | Optimal dilution of the recombinant PRL antibody should be determined by the researcher.                     |
| <b>Formula</b>             | 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide                                      |
| <b>Isotype</b>             | Mouse IgG1, kappa                                                                                            |
| <b>Clonality</b>           | Recombinant                                                                                                  |
| <b>Purity</b>              | Protein A/G affinity                                                                                         |
| <b>Uniprot ID</b>          | <b>P01236</b>                                                                                                |
| <b>Hazard Information</b>  | This recombinant PRL antibody is available for research use only.                                            |
| <b>Dilution Range</b>      | Immunohistochemistry (FFPE): 1-2ug/ml                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                              |



IHC staining of FFPE human pituitary tis...



Negative control: IHC staining of FFPE h...



SDS-PAGE analysis of purified, BSA-free ...