



Product Datasheet

LEF1 antibody (orb1151564)


Description	Rabbit monoclonal antibody to LEF1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 100-200 was used as the immunogen for the LEF1 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the LEF1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the LEF1 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	Q9UJU2
Hazard Information	This LEF1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 2-4ug/ml
Expiration Date	12 months from date of receipt.



IHC staining
of FFPE
human lymph
node wi...



IHC staining
of FFPE
human tonsil
tissue...



SDS-PAGE
analysis of
purified, BSA-
free ...