



Product Datasheet

EGF Receptor antibody (orb1151569)


Description	Mouse monoclonal antibody to EGF Receptor
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Amino acids 25-645 from the extracellular domain of the human protein were used as the immunogen for the recombinant EGF Receptor antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the recombinant EGF Receptor antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant EGF Receptor antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a, kappa
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	P00533
Hazard Information	This recombinant EGF Receptor antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human placental tis...



Negative control: IHC staining of FFPE h...



SDS-PAGE analysis of purified, BSA-free ...