



Product Datasheet

TIF1 alpha antibody (orb1151590)

Description	Mouse monoclonal antibody to TIF1 alpha
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF
Immunogen	Recombinant full-length human Transcription intermediary factor 1-alpha protein was used as the immunogen for the TIF1 alpha antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the TIF1 alpha antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the TIF1 alpha antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Protein A/G affinity
Uniprot ID	O15164
Hazard Information	This TIF1 alpha antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml
Expiration Date	12 months from date of receipt.




Immunofluorescent staining of PFA-fixed ...



FACS staining of PFA-fixed human HeLa ce...

Human Protein Microarray Specificity Validation



Analysis of HuProt (TM) microarray conta...