

Product Datasheet

Ubiquitin-conjugating enzyme E2 B antibody (orb1151595)

Description

Mouse monoclonal antibody to Ubiquitin-


Species/Host

Mouse

Reactivity

Human


Conjugation

Unconjugated

Tested Applications

FACS, IF

Immunofluorescent staining of PFA-fixed ...



FACS staining of PFA-fixed human HeLa ce...

Immunogen

Recombinant full-length human UBE2B protein was used as the immunogen for the UBE2B antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Aliquot the UBE2B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Optimal dilution of the UBE2B antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2b

Clonality

Monoclonal

Purity

Protein A/G affinity

Uniprot ID

P63146

Hazard Information

This UBE2B antibody is available for research use only.


Dilution Range

Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml

Expiration Date

12 months from date of receipt.

Human Protein Microarray Specificity Validation



Analysis of HuProt (TM) microarray conta...