



Product Datasheet

SIRT2 antibody (orb1151614)


Description	Mouse monoclonal antibody to SIRT2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF, WB
Immunogen	Recombinant full-length human SIRT2/Sirtuin 2 protein was used as the immunogen for the Sirtuin 2 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the Sirtuin 2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the Sirtuin 2 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Protein A/G affinity
Uniprot ID	Q8IXJ6
Hazard Information	This Sirtuin 2 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 1-2ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescent staining of PFA-fixed ...



Immunofluorescent staining of PFA-fixed ...



FACS staining of PFA-fixed human MCF-7 c...