


Product Datasheet


EIF4A2 antibody (orb1151650)

Description	Mouse monoclonal antibody to EIF4A2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant full-length human EIF4A2 protein was used as the immunogen for the EIF4A2 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the EIF4A2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the EIF4A2 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Protein A/G affinity
Uniprot ID	Q14240
Hazard Information	This EIF4A2 antibody is available for research use only.
Dilution Range	ELISA (order BSA-free format for coating)
Expiration Date	12 months from date of receipt.

Human Protein Microarray Specificity Validation



Analysis of HuProt (TM) microarray conta...



SDS-PAGE analysis of purified, BSA-free ...