



Product Datasheet

LMO2 antibody (orb1151693)


Description	Rabbit monoclonal antibody to LMO2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF
Immunogen	A portion of amino acids 23-140 was used as the immunogen for the recombinant Rhombotin 2 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the recombinant Rhombotin 2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant Rhombotin 2 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	P25791
Hazard Information	This recombinant Rhombotin 2 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml
Expiration Date	12 months from date of receipt.



FACS staining of K562 cells using recomb...



Immunofluorescent staining of human K562...



SDS-PAGE analysis of purified, BSA-free ...