



Product Datasheet

TdT antibody (orb1151812)


Description	Rabbit monoclonal antibody to TdT
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 1-100 from the human protein was used as the immunogen for the DNA nucleotidylexotransferase antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the DNA nucleotidylexotransferase antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the DNA nucleotidylexotransferase antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Purity	Protein A/G affinity
Uniprot ID	P04053
Hazard Information	This DNA nucleotidylexotransferase antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thymus with D...



IHC staining of FFPE human thymus tissue...



Negative control: IHC staining of FFPE h...